


Sepax UHPLC Column Products

Silica Based Packings

Stationary Phases



Description

Sepax UHPLC silica packings use high purity (<10 ppm metals), spherical silica with the particle selection of 1.8 and 2.2 μm and the pore size of 120 \AA . Their bonded phases include C18, C8, C4, Phenyl, Amino, Cyano, strong anion, strong cation, and HILIC. Their unique mono-functional bonding chemistry allows high efficiency, high selectivity and high speed separation. All sub-2 micron particles and bonded phases have excellent resistance to high pressure (>10,000 psi). As an example in Figure 1, 2.2 μm GP-C18 shows a reduced plate height of 2.22 μm , which is equivalent to plate number of 200,448 per meter.

Specifications

Phases	Particle size (μm)	Pore size (\AA)	Surface area (m^2/g)	Carbon loading	pH range
GP-C18	1.8, 2.2	120	300	17.0%	2-8.5
BR-C18	1.8, 2.2	120	300	19.0%	1.5-10.5
GP-C8	1.8, 2.2	120	300	11.0%	2-8.5
GP-C4	1.8, 2.2	120	300	8.0%	2-8.5
GP-Phenyl	1.8, 2.2	120	300	11.0%	2-8.5
HP-Cyano	1.8, 2.2	120	300	7.0%	2-8.5
HP-Amino	1.8, 2.2	120	300	4.0%	2-8.5
HP-SCX	1.8, 2.2	120	300	11.0%	2-8.5
HP-SAX	1.8, 2.2	120	300	16.0%	2-8.5
HP-Silica	1.8, 2.2	120	300	0.0%	2-8.5
HILIC Polar-100	1.8, 2.2	120	300	11.0%	2-8.5




Figure 1. Reduced plate height (h) vs. linear flow rate (v).
Column: GP-C18 (2.2 μm , 4.6x50 mm)
Mobile phase: 70% MeOH/30% H₂O
Temperature: 30 $^{\circ}\text{C}$
Detection: UV 254 nm
Test compound: naphthalene
Injection: 5 μL

Characteristics

- Most comprehensive selection of the stationary phases
- Highly controlled chemistry of monolayer formation and end-capping
- High column-to-column reproducibility
- High mechanical stability to resist the pressure as high as up to 10,000 psi
- High resolution, efficiency and selectivity for separations
- Suitable for separations of acidic, neutral and basic compounds, peptides, and proteins

Applications

- UHPLC separations
- Pharmaceuticals
- Peptides
- Proteins

Sepax sub-2 μm particles and columns can work well with both regular HPLC and UHPLC systems. Figure 2 shows an example of separation of a mixture of test compounds. The flow rate is 0.43 mL/min, generating 3060 psi for a 2.2 μm GP-C18 (3.0x150 mm) column. The flow rate was not optimized to reach minimum reduced plate height. However, it showed high separation efficiency and selectivity.




Figure 2. Separation of a mixture by regular HPLC system. Column: GP-C18 (2.2 μm , 3.0x150mm). Mobile phase: 60% Acetonitrile/40% H₂O, 0.43 mL/min. Temperature: ambient. Detection: UV 254 nm. Injection: 2 μL .

Separation of a mixture of neutral and basic compounds




Figure 3. Separation of a mixture by regular HPLC system. Column: GP-C18 (2.2 μm , 3.0x150mm). Mobile phase: 60% Acetonitrile/40% H₂O, 0.43 mL/min. Temperature: ambient. Detection: UV 254 nm. Injection: 2 μL .

Separation of barbital




Figure 4. Barbital separated by a 2.2 μm GP-C18 column and a 3 μm commercial C18 column with regular HPLC system. Columns: 4.6x50 mm. Mobile phase: 50% MeOH/50% H₂O. Flow rate: 1.0 mL/min. Temperature: 30°C. Detection: UV 214 nm. Injection: 3 μL . Sample: 1. barbital, 2. phenobarbital, 3. aprobarbital, 4. butabarbital, 5. mephobarbital, 6. pentobarbital, and 7. secobarbital.

Separation of Vitamin D Isomers




Figure 5. Vitamin D2 and D3 separated by a Sepax 2.2 µm GP-C18 column, a 2.2 µm and a 3 µm C18 columns from other vendors. Columns: 4.6x50 mm. Injection: 3 µL.

Separation of Caffeine & Metabolites






Figure 6. Caffeine and metabolites separated by a Sepax 2.2 µm GP-C18 column, a 2.2 µm and a 3 µm C18 columns from other vendors. Columns: 4.6x50mm. Mobile phase A: 0.1% NH₄OAc, pH 4.5; B: ACN. Gradient: 98%-93%-85%-85% A (0-0.75-1.50-2.25min). Flow rate: 3.0 mL/min. Temperature: 30°C. Detection: UV 275nm. Injection: 5 µL.

PolyRP Phases



Specifications

PS/DVB Particles: spherical, 80% cross-linking
 Pore size: non-porous
 Particle size: 1.0, 1.8 µm
 Surface area: <10 m²/g
 Phase structure: phenyl group
 Separation mechanism: hydrophobic interaction



Description



PolyRP phases are made of 80% cross-linking PS/DVB spherical particles. Those high rigid particles have both non-porous and porous structures with the particle selection of 1 and 1.8 µm. The phase structure is phenyl functional group that enables hydrophobic interaction. Compared with silica based reversed phases, PolyRP phases have advantages over applications at extreme pH (1-14) with special selectivity and slightly lower separation efficiency.

Characteristics

- High resolution and efficiency
- Different selectivity as the silica based reversed phases
- High column-to-column reproducibility
- Unique PolyRP reversed phases for pH 1-14
- High mechanical stability to resist the pressure as high as up to 10,000 psi
- Suitable for separations of acidic, neutral and basic compounds, peptides, and proteins

Proteomix Ion-exchange phases

Stationary phases



Specifications

Proteomix Phases	Particle size (µm)	Pore size	Functional Group	Dynamic loading capacity	pH range
SCX-NP	1.8	Non-porous	Sulfonate	65 mg/mL	2-12
WCX-NP	1.8	Non-porous	Carboxylate	25 mg/mL	2-12
SAX-NP	1.8	Non-porous	Quaternary ammonium	43 mg/mL	2-12
WAX-NP	1.8	Non-porous	Tertiary amine	35 mg/mL	2-12

Description

Proteomix ion-exchange phases have been innovatively developed for achieving high ion-exchange capacity, a breakthrough technology for non-porous resins. Proteomix packing support is composed of a rigid, spherical, highly cross-linked poly(styrene divinylbenzene) (PS/DVB) non-porous bead. The PS/DVB resin surface is grafted with a highly hydrophilic, neutral polymer thin layer with the thickness in the range of nanometers, which shields the hydrophobic PS/DVB core and eliminates their non-specific bindings with biological analytes, leading to high efficiency and high recovery separations for biological molecules. On the top of the hydrophilic layer, ion-exchange functional groups are attached via a proprietary chemistry, resulting in high capacity ion-exchange layer.

Characteristics

- Pioneered high capacity sub-2 µm non-porous ion-exchange packings
- Unprecedented high efficiency and resolution
- High column-to-column reproducibility

- High mechanical stability to resist the pressure as high as up to 10,000 psi
- High speed, high recovery and high stability

Applications

The uniqueness of non-porous Proteomix ion-exchange phases is that they combine increased capacity and resolution with their intrinsic advantages of high efficiency and high separation speed. Proteomics ion-exchange phases are especially suitable for high resolution, high efficiency and high recovery separations of proteins, oligonucleotides, peptides, polysaccharides, cell lysates, nanoparticles and nanotubes.

- UHPLC separations
- Proteins
- DNA and oligonucleotides
- Peptides
- Carbohydrates
- Cell lysates

Separation of proteins

Extremely high efficiency separation of proteins was achieved by a 1.8 μm Proteomix WCX column, as shown in Figure 1.




Figure 1. Separation of proteins by a Proteomix WCX-NP1.8 (1.8 μm , 4.6x50mm). Mobile phase: (A) 10 mM phosphate, pH 6.0; (B) A + 1.0 M NaCl; 0-75%B in 25 min; Flow rate: 0.75 mL/min. Temperature: ambient. Detection: UV 280 nm. Sample: (1) Ribonuclease A, (2) Cytochrome C and (3) Lysozyme. Injection: 5 μL .

Separation of ovalbumin major components from its impurities by anionic exchange column

Ovalbumin from Sigma has two major components and a number of minor components or impurities. Resolution of those components is a good indication of the efficiency and resolving power of a column. A short Proteomix SAX-NP1.8 column (1.8 μm , 4.6x30mm) very well resolved those components, as shown in Figure 2.




Figure 2. Separation of multiple components of ovalbumin by a Proteomix SAX-NP1.8 (1.8 μm , 4.6x30mm). Mobile phases: (A) 20 mM Tris buffer, pH 8.0; (B) A + 0.5 M NaCl. Gradient: 0-75%-100%B (0-21-21.5min). Flow rate: 0.4 mL/min. Temperature: ambient. Backpressure: 1,800 psi. Detection: UV 214 nm. Injection: 5 μL (10 mg/mL of ovalbumin in H₂O).

Ordering Information

1.8 μm Reversed and Normal Phases (ID x length mm)						
Phases	2.1x30	2.1x50	2.1x100	2.1x150	3.0x30	3.0x50
GP-C18 (120 \AA)	101181-2103	101181-2105	101181-2110	101181-2115	101181-3003	101181-3005
BR-C18 (120 \AA)	102181-2103	102181-2105	102181-2110	102181-2115	102181-3003	102181-3005
GP-C8 (120 \AA)	107081-2103	107081-2105	107081-2110	107081-2115	107081-3003	107081-3005
GP-C4 (120 \AA)	109041-2103	109041-2105	109041-2110	109041-2115	109041-3003	109041-3005
GP-Phenyl (120 \AA)	111361-2103	111361-2105	111361-2110	111361-2115	111361-3003	111361-3005
HP-CN (120 \AA)	113311-2103	113311-2105	113311-2110	113311-2115	113311-3003	113311-3005
HP-NH ₂ (120 \AA)	115301-2103	115301-2105	115301-2110	115301-2115	115301-3003	115301-3005
HP-SCX (120 \AA)	120361-2103	120361-2105	120361-2110	120361-2115	120361-3003	120361-3005
HP-Silica (120 \AA)	117001-2103	117001-2105	117001-2110	117001-2115	117001-3003	117001-3005
HILIC Polar-100	131581-2103	131581-2105	131581-2110	131581-2115	131581-3003	131581-3005

1.8 μm Reversed and Normal Phases (ID x length mm)						
Phases	3.0x100	3.0x150	4.6x30	4.6x50	4.6x100	4.6x150
GP-C18 (120 Å)	101181-3010	101181-3015	101181-4603	101181-4605	101181-4610	101181-4615
GP-C18 (120 Å)	102181-3010	102181-3015	102181-4603	102181-4605	102181-4610	102181-4615
GP-C8 (120 Å)	107081-3010	107081-3015	107081-4603	107081-4605	107081-4610	107081-4615
GP-C4 (120 Å)	109041-3010	109041-3015	109041-4603	109041-4605	109041-4610	109041-4615
GP-Phenyl (120 Å)	111361-3010	111361-3015	111361-4603	111361-4605	111361-4610	111361-4615
HP-CN (120 Å)	113311-3010	113311-3015	113311-4603	113311-4605	113311-4610	113311-4615
HP-NH ₂ (120 Å)	115301-3010	115301-3015	115301-4603	115301-4605	115301-4610	115301-4615
HP-SCX (120 Å)	120361-3010	120361-3015	120361-4603	120361-4605	120361-4610	120361-4615
HP-Silica (120 Å)	117001-3010	117001-3015	117001-4603	117001-4605	117001-4610	117001-4615
HILIC Polar-100	131581-3010	131581-3015	131581-4603	131581-4605	131581-4610	131581-4615

2.2 μm Reversed and Normal Phases (ID x length mm)						
Phases	2.1x30	2.1x50	2.1x100	2.1x150	3.0x30	3.0x50
GP-C18 (120 Å)	101182-2103	101182-2105	101182-2110	101182-2115	101182-3003	101182-3005
GP-C18 (120 Å)	102182-2103	102182-2105	102182-2110	102182-2115	102182-3003	102182-3005
GP-C8 (120 Å)	107082-2103	107082-2105	107082-2110	107082-2115	107082-3003	107082-3005
GP-C4 (120 Å)	109042-2103	109042-2105	109042-2110	109042-2115	109042-3003	109042-3005
GP-Phenyl (120 Å)	111362-2103	111362-2105	111362-2110	111362-2115	111362-3003	111362-3005
HP-CN (120 Å)	113312-2103	113312-2105	113312-2110	113312-2115	113312-3003	113312-3005
HP-NH ₂ (120 Å)	115302-2103	115302-2105	115302-2110	115302-2115	115302-3003	115302-3005
HP-SCX (120 Å)	120362-2103	120362-2105	120362-2110	120362-2115	120362-3003	120362-3005
HP-Silica (120 Å)	117002-2103	117002-2105	117002-2110	117002-2115	117002-3003	117002-3005
HILIC Polar-100	131582-2103	131582-2105	131582-2110	131582-2115	131582-3003	131582-3005
Phases	3.0x100	3.0x150	4.6x30	4.6x50	4.6x100	4.6x150
GP-C18 (120 Å)	101182-3010	101182-3015	101182-4603	101182-4605	101182-4610	101182-4615
BR-C18 (120 Å)	102182-3010	102182-3015	102182-4603	102182-4605	102182-4610	102182-4615
GP-C8 (120 Å)	107082-3010	107082-3015	107082-4603	107082-4605	107082-4610	107082-4615
GP-C4 (120 Å)	109042-3010	109042-3015	109042-4603	109042-4605	109042-4610	109042-4615
GP-Phenyl (120 Å)	111362-3010	111362-3015	111362-4603	111362-4605	111362-4610	111362-4615
HP-CN (120 Å)	113312-3010	113312-3015	113312-4603	113312-4605	113312-4610	113312-4615
HP-NH ₂ (120 Å)	115302-3010	115302-3015	115302-4603	115302-4605	115302-4610	115302-4615
HP-SCX (120 Å)	120362-3010	120362-3015	120362-4603	120362-4605	120362-4610	120362-4615
HP-Silica (120 Å)	117002-3010	117002-3015	117002-4603	117002-4605	117002-4610	117002-4615
HILIC Polar-100	131582-3010	131582-3015	131582-4603	131582-4605	131582-4610	131582-4615

1.8 μm PolyRP-NP phase (ID x length mm)						
Phases	2.1x30	2.1x50	2.1x100	2.1x150	3.0x30	3.0x50
PolyRP-NP1.8	262002-2103	262001-2105	262001-2110	262001-2115	262001-3003	262001-3005
Phases	3.0x100	3.0x150	4.6x30	4.6x50	4.6x100	4.6x150
PolyRP-NP1.8	262002-3010	262001-3015	262001-4603	262001-4605	262001-4610	262001-4615

1.8 μm Proteomix Ion-exchange phases (ID x length mm)						
Phases	2.1x30	2.1x50	2.1x100	2.1x150	3.0x30	3.0x50
SCX-NP1.8	401NP2-2103	401NP2-2105	401NP2-2110	401NP2-2115	401NP2-3003	401NP2-3005
WCX-NP1.8	402NP2-2103	402NP2-2105	402NP2-2110	402NP2-2115	402NP2-3003	402NP2-3005
SAX-NP1.8	403NP2-2103	403NP2-2105	403NP2-2110	403NP2-2115	403NP2-3003	403NP2-3005
WCX-NP1.8	404NP2-2103	404NP2-2105	404NP2-2110	404NP2-2115	404NP2-3003	404NP2-3005
Phases	3.0x100	3.0x150	4.6x30	4.6x50	4.6x100	4.6x150
SCX-NP1.8	401NP2-3010	401NP2-3015	401NP2-4603	401NP2-4605	401NP2-4610	401NP2-4615
WCX-NP1.8	402NP2-3010	402NP2-3015	402NP2-4603	402NP2-4605	402NP2-4610	402NP2-4615
SAX-NP1.8	403NP2-3010	403NP2-3015	403NP2-4603	403NP2-4605	403NP2-4610	403NP2-4615
WCX-NP1.8	404NP2-3010	404NP2-3015	404NP2-4603	404NP2-4605	404NP2-4610	404NP2-4615