

Sepax SRT-10C SEC Media for High Throughput Preparative Size Exclusion Chromatography

Size Exclusion Chromatography (SEC) is a chromatographic technique that allows for separation of molecules based on their hydrodynamic size. In addition to being a common method for analytical characterization and quality control analysis of biological drug candidates, SEC is also a technique well employed in downstream biomolecule preparative purification (more commonly aggregate removal). Sepax Technologies, Inc. is a leading manufacturer of liquid chromatography (LC) related resins, consumables, and columns. The Sepax SRT-10C SEC column for Fast Purification, offers a robust, hydrophilic surface modified rigid silica matrix that allows for faster flow rates translating to shorter run times.

In the following paper, *Development of a robust and semi-automated two-step antibody purification process*, recently published by Eli Lilly, Yang et al. highlight the advantages of the Sepax SRT-10C SEC column chemistry and propose a quick but reliable cleaning method for Endotoxin removal.

Literature reference: Yang, Xiaomin, et al. "Development of a robust and semi-automated two-step antibody purification process." *Mabs*. Vol. 13. No. 1. Taylor & Francis, 2021.


Table 1. Performance comparison of SRT-10 C SEC-300 compared to superdex 200 column.

| Parameter / Results Column | Loading Vol. (mL) | Flow Rate (mL/min) | Running Time (min) | Resolution* | Effluent Vol (mL) | Effluent Vol Increase | Yield % | aSEC Purity % |
|---------------------------------------|-------------------|--------------------|--------------------|-------------|-------------------|-----------------------|---------|---------------|
| Superdex 200 | 6 | 1.5 | 80 | 1.05 | 16.0 | 2.7-fold | 76.3 | 97.1 |
| (16 x 600mm, 120 mL) | 6 | 1.5 | 80 | 1.05 | 16.0 | 2.7-fold | 75.9 | 97.1 |
| SRT-10 C SEC-300 (30 x 300mm, 210 mL) | 9 | 7.5 | 40 | 1.08 | 25.5 | 2.8-fold | 79.1 | 97.1 |
| | 9 | 7.5 | 40 | 1.07 | 25.5 | 2.8-fold | 79.5 | 97.0 |

*Resolution was calculated from $[(\text{Volume}_{\text{peak2}} - \text{Volume}_{\text{peak1}}) / ((\text{Width}_{1/2, \text{peak2}} + \text{Width}_{1/2, \text{peak1}}) / 2)]$ by UNICORN software in AKTA pure.

Table 1, Sepax SRT-10C SEC-300 column showed the advantages of faster flow rate, shorter running time without sacrificing the resolution, yield and % purity of the antibody sample.

Figure 1 below illustrates that high resolution is achieved as per the red trace below. Based on the conditions outlined below, it is conceivable that 1-2% CV loading with a sample range of 2-7 mg/mL at 5 mL/min could be used to maximize loading while not drastically effecting resolution (or provide sufficient purification based on a acceptable range) .The blue trace below represents extreme conditions with a 25X load at nearly double the flow rate.



From Table 2, endotoxin contamination was significantly reduced to less than 0.5 EU/mL in washes which IPA concentration was greater than 25%. Prep SEC column can be washed and stored in 25% IPA during routine operation and completely cleaned with 70% IPA in instances of endotoxin contamination.

Table 2. Endotoxin reduction following clean-in-place using various concentrations of isopropyl alcohol (IPA).

| Recovery % IPA% | Pre-IPA Protein Recovery% | Pre-IPA Endotoxin Recovery% | Post-IPA Protein Recovery% | Post-IPA IPA Total EU | Endotoxin Reduction % |
|--------------------|---------------------------|-----------------------------|----------------------------|-----------------------|-----------------------|
| 0 | 95 | 107 | 98 | 220* | N/A |
| 10 | 101 | 131 | 98 | 121 | 45 |
| 25 | 99 | 83 | 103 | 11 | 95 |
| 50 | 99 | 97 | 100 | 2 | 99 |
| 70 | 97 | 102 | 100 | 5 | 98 |

*baseline endotoxin level, which was used to calculate the endo removal%.

Endotoxin Wash: To establish an effective cleaning procedure for the SRT-10C column with IPA aqueous solution, a purified control antibody at 2 mg/mL alone and a mixture of antibody with bovine thyroglobulin at 3 mg/ml containing 3137 EU/mL endotoxin were used for testing. After the antibody/thyroglobulin mixture passes through the column, it was first washed with one CV of 70% IPA solution followed by one CV of PBS before injecting the pure antibody sample. The eluents from both injections were collected for concentration and endotoxin measurement; the protein and endotoxin recovery were calculated relative to the loading amount. This experiment

Order Information:

SRT-10C SEC-300, 10 μm , 300 \AA , 30 x 300 mm
Part Number: [239300-30030](#)

See our full line up of Sepax SRT-10C SEC Columns [here](#).
Ranging from 100 \AA to 2000 \AA in Pore Sizes and Column Volumes up to ~590 mL

Sepax SRT-10C SEC-100
Sepax SRT-10C SEC-300
Sepax SRT-10C SEC-500
Sepax SRT-10C SEC-1000
Sepax SRT-10C SEC-2000

Request a [quote](#) today!

**SPECIALIZED & INNOVATIVE
CHROMATOGRAPHY EXPERTS**



Visit our website: Sepax-Tech.com
Contact us: 1-877-SEPAX-US and Info@sepax-tech.com
Reach our Technical Support: TechSupport@sepax-tech.com

FOLLOW US

