

SP1022- Fast Protein Purification with 5 cm ID Large Preparative Size SRT-10 SEC Column

HPLC Column, Sepax, SRT-10 SEC-300, 10um, 300 A 50 x 250mm
(225300-50025)



Protein purification with SRT10 SEC-300 50 x 250 mm

Column: SRT10 SEC-300 (10 μm , 300 \AA , 50 x 250 mm) **Detection:** 280 nm
Mobile phase: 150 mM sodium phosphate buffer, pH 7.0 **Column temperature:** Ambient
Flow rate: 25 mL/min **Injection Volume:** 10 mL 10 mg/mL BSA in water
Instrument: AKTA Pure

